



\* Cross Reactivity  
Cross reactivity of SARS-CoV-2/FluA/FluB+ADV/RSV/SA Antigen Combo Rapid Test Kit (LFIA) was evaluated by testing commercial and pathogenic microorganisms listed in the following table that may be present in the clinical samples. Each of the bacterium, viruses, and yeast were tested in triplicate with no false positive results of SARS-CoV-2 virus, Influenza A, Influenza B, ADV, RSV and STREP A.

Potential Cross-Reactant	Concentration Tested	Cross-Reactivity (Yes/No)					
		SARS-CoV-2	Flu A	Flu B	ADV	RSV	STREP A
Human coronavirus OC43	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	No	No	No	No
Human coronavirus NL33	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	No	No	No	No
Human coronavirus HKU1	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	No	No	No	No
Human coronavirus 229E	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	No	No	No	No
MERS-coronavirus	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	No	No	No	No
SARS-coronavirus	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	No	No	No	No
SARS-CoV-2	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	/	No	No	No	No	No
H1N1(2009)	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	/	No	No	No	No
Influenza A H1N1 Seasonal	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	/	No	No	No	No
Influenza A H3N2	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	/	No	No	No	No
Influenza A H7N9	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	/	No	No	No	No
Influenza B Victoria	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	/	No	No	No	No
Influenza B Yamagata	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	/	No	No	No	No
Parainfluenza virus Type 1	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	No	No	No	No
Enterovirus CA16e	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	No	No	No	No
Adenovirus	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	/	No	No	No
ADV-2	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	/	No	No	No
ADV-3	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	/	No	No	No
ADV-4	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	/	No	No	No
ADV-5	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	/	No	No	No
RSV-A	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	No	No	/	No
RSV-B	1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL	No	No	No	No	/	No
<i>Mycoplasma pneumoniae</i>	1.0 x 10 <sup>3</sup> CFU/mL	No	No	No	No	No	No
<i>Staphylococcus aureus</i>	1.0 x 10 <sup>3</sup> CFU/mL	No	No	No	No	No	No
<i>Staphylococcus epidermidis</i>	1.0 x 10 <sup>3</sup> CFU/mL	No	No	No	No	No	No
<i>Bordetella pertussis</i>	1.0 x 10 <sup>3</sup> CFU/mL	No	No	No	No	No	No
<i>Legionella pneumophila</i>	1.0 x 10 <sup>3</sup> CFU/mL	No	No	No	No	/	No
<i>Streptococcus pneumoniae</i>	1.0 x 10 <sup>3</sup> CFU/mL	No	No	No	No	No	No
<i>Haemophilus influenzae</i>	1.0 x 10 <sup>3</sup> CFU/mL	No	No	No	No	No	No
<i>Mycobacterium tuberculosis</i>	1.0 x 10 <sup>3</sup> CFU/mL	No	No	No	No	No	No
<i>Candida albicans</i>	1.0 x 10 <sup>3</sup> CFU/mL	No	No	No	No	No	No

\* Interfering Substances Effect  
A study was performed to evaluate and demonstrate that the endogenous substances naturally present or drugs that may be artificially introduced into clinical samples do not inference with the detection of SARS-CoV-2, Influenza A, Influenza B, ADV, RSV and STREP A in the SARS-CoV-2/Flu A/H1N1/RSV/SA Antigen Combo Rapid Test Kit (LFIA) at the concentrations listed below. Dilute test items with Swab as sample matrix in the absence or presence of heat inactivated SARS-CoV-2, Influenza A, Influenza B, ADV, RSV and STREP A.

Type	Potential Interfering Substances	Concentration	Interference(Yes/No)
Endogenous Substance	Mucin	2% w/v	No
	Whole Blood	5% w/v	No
	Icteric (Bilirubin)	40 mg/dL	No
	Rheumatoid factor	200 IU/mL	No
	Triglycerides	1.5 mg/L	No
	Hemoglobin	100 mg/L	No
	Anti-nuclear antibody	>1:40	No
	Total IgG	90 g/L	No
	Total IgM	4 g/L	No
	Total IgA	80 g/L	No
Exogenous Substance	Tamiflu (Osetamivir Phosphate)	0.55% w/v	No
	Fluticasone Propionate	5% w/v	No
	Flicozazole	5% w/v	No
	Zincum gluconium (I.e., Zicam)	5% w/v	No
	Alkalol	10% w/v	No
	Phenylephrine hydrochloride	15% v/v	No
	Oxymentazolin Hydrochloride	15% v/v	No
	Cromolyn	15% w/v	No
	Oxymetazoline	15% w/v	No
	Galiphimia glauca, Sabadilla	20% w/v	No
Exogenous Substance	Albuterol	0.005 mg/dL	No
	Acarbose	0.03 mg/dL	No
	Osetamivir	0.04 mg/dL	No
	Chlorpheniramine	0.08 mg/dL	No
	Diphenhydramine	0.08 mg/dL	No
	Glimeprid (Sulfonylureas)	0.164 mg/dL	No
	Chlorothiazide	2.7 mg/dL	No
	Acetylsalicylic acid	3 mg/dL	No
	Amoxicillin	5.4 mg/dL	No
	Ibuprofen	21.9 mg/dL	No
Exogenous Substance	Beclometason	4.79 mg/mL	No
	Indapamide	140 ng/mL	No
	Flunisolide	0.61 µg/mL	No
	Guaiacol glyceryl ether	120 µg/mL	No
	Zanamivir	17.3 µg/mL	No
	Tobramycin	24.03 µg/mL	No
	Sulfur	9.23 µg/mL	No
	Ribavirin	26.7 µg/mL	No
	Ephedrine	0.1 mg/mL	No
	Meprobamate	0.13 mg/mL	No
Exogenous Substance	Menthol	0.15 mg/mL	No
	Budesonide	0.5 mg/mL	No
	Triamcinolone	0.8 mg/mL	No
	Dexamethasone	0.8 mg/mL	No
	Lopinavir	16.4 µg/L	No
	Ritonavir	16.4 µg/L	No
	Chloroquine phosphate	0.99 mg/L	No
	Ivermectin	4.4 mg/L	No

**Clinical performance**  
\* Clinical performance with throat swabs  
The performance of SARS-CoV-2/FluA/FluB+ADV/RSV/SA Antigen Combo Rapid Test Kit (LFIA) was established with 134 throat swabs collected from patients with COVID-19 symptoms within 7 days after onset of symptoms. Two swabs were collected with the same people, an throat swab tested directly using SARS-CoV-2/FluA/FluB+ADV/RSV/SA Antigen Combo Rapid Test Kit (LFIA) and a throat swab tested by the RT-PCR Test Kit. Clinical samples were evaluated to be positive or negative using RT-PCR reference method.

**2. Influenza A Test**  
The performance of SARS-CoV-2/FluA/FluB+ADV/RSV/SA Antigen Combo Rapid Test Kit (LFIA) was established with 103 throat swabs collected from patients with Influenza symptoms within 7 days after onset of symptoms. Two swabs were collected with the same people, an throat swab tested directly using SARS-CoV-2/FluA/FluB+ADV/RSV/SA Antigen Combo Rapid Test Kit (LFIA) and a throat swab tested by the Comparison Test Kit. Clinical samples were evaluated to be positive or negative using Comparison reference method.

**3. Influenza B Test**  
The performance of SARS-CoV-2/FluA/FluB+ADV/RSV/SA Antigen Combo Rapid Test Kit (LFIA) was established with 106 throat swabs collected from patients with ADV symptoms within 7 days after onset of symptoms. Two swabs were collected with the same people, an throat swab tested directly using SARS-CoV-2/FluA/FluB+ADV/RSV/SA Antigen Combo Rapid Test Kit (LFIA) and a throat swab tested by the Comparison Test Kit. Clinical samples were evaluated to be positive or negative using Comparison reference method.

**4. ADV Test**  
The performance of SARS-CoV-2/FluA/FluB+ADV/RSV/SA Antigen Combo Rapid Test Kit (LFIA) was established with 121 throat swabs collected from patients with ADV symptoms within 14 days after onset of symptoms. Two swabs were collected with the same people, an throat swab tested directly using SARS-CoV-2/FluA/FluB+ADV/RSV/SA Antigen Combo Rapid Test Kit (LFIA) and a throat swab tested by the Comparison Test Kit. Clinical samples were evaluated to be positive or negative using Comparison reference method.

**5. RSV Test**  
The performance of SARS-CoV-2/FluA/FluB+ADV/RSV/SA Antigen Combo Rapid Test Kit (LFIA) was established with 146 throat swabs collected from patients with ADV symptoms within 14 days after onset of symptoms. Two swabs were collected with the same people, an throat swab tested directly using SARS-CoV-2/FluA/FluB+ADV/RSV/SA Antigen Combo Rapid Test Kit (LFIA) and a throat swab tested by the Comparison Test Kit. Clinical samples were evaluated to be positive or negative using Comparison reference method.

**6. STREP A Test**  
The performance of SARS-CoV-2/FluA/FluB+ADV/RSV/SA Antigen Combo Rapid Test Kit (LFIA) was established with 118 throat swabs collected from patients with ADV symptoms within 14 days after onset of symptoms. Two swabs were collected with the same people, an throat swab tested directly using SARS-CoV-2/FluA/FluB+ADV/RSV/SA Antigen Combo Rapid Test Kit (LFIA) and a throat swab tested by the Comparison Test Kit. Clinical samples were evaluated to be positive or negative using Comparison reference method.

**Warnings and Precautions**  
\* This test kit is used for in vitro diagnosis only.  
\* This test kit should be used by qualified person with professional experience or proper training.  
\* This test kit should not be used if the package and samples for transport media will reduce sensitivity. The test cassette should not be used if being wet or polluted.  
\* Proper protection should be taken during testing to avoid splashing when adding sample.  
\* Dispose of all used or damaged test cassettes, sampling tubes, droppers, swabs or other kit components as biohazardous materials.  
\* Negative results do not rule out SARS-CoV-2, Influenza A, Influenza B, ADV, RSV and STREP A infection, particularly in those who have been in contact with the virus.  
\* The sample should not be thick, as thick samples may result in abnormal detection.

**References**  
1. LV Zhang, PW Chen, GW Zheng, et al. Research progress on COVID-19 in recent methods. Modern Medicine and Clinic, 2020, 35(3): 411-416.  
2. K Tugba, W Ralph, L Hahn. Molecular and Immunological Diagnostic Tests of COVID-19: Current Status and Challenges. iScience, 2020, 23 (8): DOI: 10.1016/j.isci.2020.101406  
3. WHO recommendations on the use of rapid testing for influenza diagnosis. World Health Organization, July 2005.  
4. Jian X, Fan Y, Wang C, et al. Seroprevalence of Neutralizing Antibodies against Six Human Adenovirus Types Indicates the Low Level of Herd Immunity in Young Children from Guangzhou, China[J]. Chinese Journal of Virology: English Edition, 2021, 36(3).  
5. Jang J, Chen Y, Wang Z. Specific IgG antibodies against F and G glycoproteins of respiratory syncytial virus (RSV) in asthmatic children after infection with the virus[J]. Chinex Journal of Pediatrics, 1998.  
6. Kim S, Lee NY. Asymptomatic infection by streptococcus pyogenes in schoolchildren and diagnostic of antidioxyribonuclease B. J Korean Med Sci, 2005, 20: 938~940.

## SARS-CoV-2/FluA/FluB+ADV /RSV/SA -Antigen-Kombination-Schnelltestkit (LFIA)

### Test-Kassette: 1St./Beutel

REF	Spezifikation	REF	Spezifikation
123015-01-01	1 St./Karton	123015-20-01	20 St./Karton
123015-02-01	2 St./Karton	123015-50-01	50 St./Karton
123015-05-01	5 St./Karton	123015-100-01	100 St./Karton

\* Kreuzreaktivität  
Die Kreuzreaktivität des SARS-CoV-2/FluA/FluB+ADV/RSV/SA-Antigen-Kombination-Schnelltestkits (LFIA) wurde durch das Testen von kommerziellen und pathogenen Mikroorganismen bewertet, die in der folgenden Tabelle aufgeführt sind und in den klinischen Proben vorhanden sein können. Jedes Bakterium, jeder Virus und jede Hefe wurde in dreifacher Ausführung getestet, wobei keine falsch-positive Ergebnisse für SARS-CoV-2-Virus, Influenza A, Influenza B, ADV, RSV und STREP A festgestellt wurden.

<